



SAUFLO^N
*Multi*TM

*The world's most simple to use
preservative free genuine one step multifunction solution*

Introducing

SAUFLON *Multi*[™]

*"The world's most simple
to use Preservative Free
Genuine One Step
Multifunction Solution"*



THE WORLD'S FIRST
Preservative Free Multifunction Solution.



THE WORLD'S FIRST
Multifunction Solution Based On Peroxide



THE WORLD'S FIRST
Preservative Free Multifunction Solution That
Eliminates The Need For Digital Cleaning.



THE WORLD'S FIRST
Preservative Free Genuine One Step
Multifunction Solution.



THE WORLD'S MOST
Simple To Use Preservative Free
Multifunction Solution.


Sauflon Multi™ is the solution every contact lens practitioner has been looking for - efficient, simple, safe, economical and virtually impossible for the patient to be non compliant.

Sauflon Multi™ is the solution every contact lens patient has been looking for - put the lenses in the solution at night, take them out and straight into the eyes in the morning. No cleaning, rinsing or neutralising - a genuine one step product.

Sauflon Multi™ Functions

Cleans

Sauflon Multi™ is the first multi function solution to dispense with the need for digital cleaning. The **Sauflon Multi™** formulation contains a non ionic surfactant cleaner which is activated during the early stages of neutralisation at which time a vigorous vortex cleaning of the lenses takes place.



Disinfects

The antimicrobial efficacy of 3% Hydrogen Peroxide is widely considered to be the most effective disinfectant used in contact lens systems providing highly effective disinfection against both bacteria and fungi.

Neutralises

The inclusion of a built in Platinum Disk in the lens case ensures that the Hydrogen Peroxide is completely neutralised in 6 hours or overnight.

Lubricates

Sauflon Multi™ formulation contains a lubricating agent which enhances wearer comfort.

Stores

Sauflon Multi™ is preservative free, thus eliminating the risk of preservative uptake during storage.

Removes Protein

Sauflon Multi™ also removes protein and other deposits from the lens surface.

SAUFLON *Multi*[™]

Sauflon Multi[™] is immediately available within the Sauflon *Protector Plan*[™] where three months supply of Sauflon lenses and solutions are placed in a box and designated by patient. Every three months the package either arrives at your practice for the patient to collect (*Sauflon System*) or is sent directly to your patient's home on your behalf (*Optician Direct*[™]).



Optician Direct[™]
Directly to your patient



**Sauflon System
Practice Collection**

 **SAUFLON
Protector Plan[™]**

SAUFIONTM *Multi*

PRODUCT PRESENTATION



BENEFITS TO YOU.

- *Optician Only*
- *Preservative Free*
- *Outstanding Simplicity*
- *Highly Effective*
- *Extremely Safe*
- *Ensures Patient Compliance*
- *Genuine One Step Solution*
- *Flat Bottle*
- *Competitively Priced*
- *Immediately Available In Saufion Protector PlanTM*

HEAD OFFICE

**SAUFLON PHARMACEUTICALS LIMITED,
49-53 YORK STREET, TWICKENHAM,
MIDDLESEX, TW1 3LP, ENGLAND.**

TELEPHONE: 00 44 (0)208 322 4200 FAX: 00 44 (0)208 891 3001

OTHER SAUFLON COMPANIES:

SAUFLON ITALIA SRL Telephone: 00 3906 51 91 504 Fax: 00 3906 50 41 317

SAUFLON FRANCE Telephone: 00 33 1 41 11 19 75 Fax: 00 33 1 41 11 19 75

SAUFLON MIDDLE EAST Telephone: 00 962 6 56 99 83 2 Fax: 00 962 6 56 20 67 6

SAUFLON INTERNATIONAL Telephone: 00 3906 51 95 57 34 Fax: 00 3906 50 41 317

SAUFLON RUSSIA Telephone: 00 44 (0)208 501 3066 Fax: 00 44 (0)208 501 3009

SAUFLON CENTRAL EUROPE Telephone: 00 31 16 83 56 56 5 Fax: 00 31 16 83 56 56 6

SAUFLON AUSTRALIA (PTY) LTD Telephone: 00 61 8 82 33 6312 Fax: 00 61 8 84 10 0348

SAUFLON BALTICS Telephone: 00 371 7 37 51 92 Fax: 00 371 7 37 09 87

SAUFLON GERMANY Telephone: 00 49 180 324 7346 Fax: 00 49 180 324 7346

SAUFLON PHARMACEUTICALS INC., USA Telephone: 00 1 800 635 6967 Fax: 00 1 800 644 3622



SAUFLON **Only for opticians**
The company you can trust